

## Less is more

**FullLife**—a minimal-contact, full-face mask—has less bulk while providing more comfort and convenience.

A revolutionary design makes FullLife a great choice for men and women who need a full-face mask but prefer something more compact with a clear line of sight.

We took a minimalist approach in developing FullLife. We made it smaller and lighter, yet ensured that it is durable. The conventional forehead pad was eliminated to provide a clear line of sight in the process.

We reduced leaks with a dual-textured, grooved cushion that creates a robust seal even with patient movement. Our straightforward, step-by-step numbered headgear system reduces application and removal effort while improving fit. With less noise, FullLife is more comfortable to live with. And with just one faceplate to fit all sizes – plus our FitPack and DuoPack options – FullLife is more comfortable for you as well. Helping you fit 100% of your patients, 100% of the time sometimes means giving you less ... and much, much more.

**PHILIPS**

RESPIRONICS



Elimination of the forehead pad creates a clear line of vision.

Integrated exhalation port makes operation exceptionally quiet.

One faceplate fits all 3 cushion sizes (S, M, L).

Safety entrainment valve provides quick access to room air if pressure is less than 3cm H<sub>2</sub>O.

Swivel quick-release tabs allow the hose to be easily disconnected.

Two 360° swivels allow tubing to move freely.

Branding and size indicators provide easy identification for patients and clinicians.


Chin support band ensures a secure fit and mask stability.



Step-by-step numbered headgear system makes set-up and use easy.

Dual-textured cushion surface aids assembly while preventing mask movement on the face.

Cushion "ribs" add extra stability to the cushion and help prevent leaks into the eyes.  
 Grooved cushion creates a robust seal that accommodates movement.  
 SST flap and textured cushion design enhance comfort.



Item		DuoPack (extra cushion of the same size)	
FullLife mask	Part number	Part number	
Small mask with headgear	1047939	1052167	
Medium mask with headgear	1047940	1052168	
Large mask with headgear	1047941	1052169	
FitPack (1 faceplate and 1 cushion of each size)			
FitPack, size S,M,L	1047942		
FullLife mask without exhalation port			
Small mask with headgear	1047943		
Medium mask with headgear	1047944		
Large mask with headgear	1047995		
Replacement cushions		Small	Medium
FullLife cushion	1047922	1047923	1047924
FullLife cushion without exhalation port	1047938	1047996	1048000
Parts and accessories			Part number
FullLife headgear system			1047975

Respironics, DuoPack, FitPack, and FullLife are trademarks of Respironics, Inc. and its affiliates.



©2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Broudy KB 2/1/10 MCI 4103129 PN 1065543

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific  
 +65 6882 5282  
 Respironics Australia  
 +61 (2) 9666 4444  
 Respironics Europe, Middle East, Africa  
 +33 1 47 52 30 00  
 Respironics United Kingdom  
 +44 800 1300 845

Philips Respironics  
 1010 Murry Ridge Lane  
 Murrysville, PA 15668

Customer Service  
 +1 724 387 4000  
 800 345 6443 (toll free, US only)  
 www.philips.com/respironics